

EPM2AIP1 Antibody

PACO37046

Description

This EPM2AIP1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO37046
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	EPM2AIP 1 antibody, EPM2A (laforin) interacting protein 1 antibody, EPM2A (laforin) IP1 antibody, EPM2A-interacting protein 1 antibody, Epm2aip1 antibody, EPMIP_HUMAN antibody, FLJ11207 antibody, KIAA0766 antibody, Laforin interacting protein 1 antibody, Laforin-interacting protein antibody, My007 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q7L775
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human EPM2A-interacting protein 1 protein (251-607AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human kidney tissue using PACO37046 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human heart tissue using PACO37046 at dilution of 1:100