

HSD11B1L Antibody

PACO37334

Description

This HSD11B1L Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO37334
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	HSD11B1L antibody, HSD3 antibody, SCDR10 antibody, SDR26C2Hydroxysteroid 11-beta-dehydrogenase 1-like protein antibody, EC 1.1.1.- antibody, 11-beta-hydroxysteroid dehydrogenase type 3 antibody, 11-DH3 antibody, 11-beta-HSD3 antibody, Short chain dehydrogenase/reductase family 26C member 2 antibody, Short-chain dehydrogenase/reductase 10 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q7Z5J1
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Hydroxysteroid 11-beta-dehydrogenase 1-like protein (17-87AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human prostate cancer using PACO37334 at dilution of 1:100