

## DDX42 Antibody

PACO37422

### Description

---

This DDX42 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |                                                                                                                                                                                                                                                                                                                         |
|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>          | PACO37422                                                                                                                                                                                                                                                                                                               |
| <b>Contents:</b>     | 50µg<br>Bradford Reagent: 1 vial (2ml)                                                                                                                                                                                                                                                                                  |
| <b>Category:</b>     | -                                                                                                                                                                                                                                                                                                                       |
| <b>Synonyms:</b>     | DDX42ATP-dependent RNA helicase DDX42 antibody, EC 3.6.4.13 antibody, DEAD box protein 42 antibody, RNA helicase-like protein antibody, RHELP antibody, RNA helicase-related protein antibody, RNAHP antibody, SF3b DEAD box protein antibody, Splicing factor 3B-associated 125 kDa protein antibody, SF3b125 antibody |
| <b>Clone:</b>        | Polyclonal                                                                                                                                                                                                                                                                                                              |
| <b>Applications:</b> | <b>ELISA</b> <b>IHC</b> <b>IF</b>                                                                                                                                                                                                                                                                                       |
| <b>Conjugation:</b>  | Non-conjugated                                                                                                                                                                                                                                                                                                          |
| <b>Reactivity:</b>   | Human                                                                                                                                                                                                                                                                                                                   |

### Antibody Data

---

|                           |                                                                               |
|---------------------------|-------------------------------------------------------------------------------|
| <b>Isotype:</b>           | IgG                                                                           |
| <b>Uniprot:</b>           | Q86XP3                                                                        |
| <b>Host Species:</b>      | Rabbit                                                                        |
| <b>Purification:</b>      | >95%, Protein G purified                                                      |
| <b>Immunogen:</b>         | Recombinant Human ATP-dependent RNA helicase DDX42 protein (519-819AA)        |
| <b>Immunogen Species:</b> | Homo sapiens (Human)                                                          |
| <b>Buffer:</b>            | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| <b>Form:</b>              | Liquid                                                                        |

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.



## Preparation & Storage


**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

| Recommended Dilutions: | Application | Recommended Dilution |
|------------------------|-------------|----------------------|
|                        | IHC         | 1:20-1:200           |
|                        | IF          | 1:50-1:500           |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

### Image



### Description

Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO37422 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human lung tissue using PACO37422 at dilution of 1:100

Immunofluorescence staining of HeLa cells with PACO37422 at 1:200, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum.

The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).