

ARMC9 Antibody

PACO37430

Description

This ARMC9 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO37430
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	ARMC9 antibody, KIAA1868LisH domain-containing protein ARMC9 antibody, Armadillo repeat-containing protein 9 antibody, Melanoma/melanocyte-specific tumor antigen KU-MEL-1 antibody, NS21 antibody
Clone:	Polyclonal
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q7Z3E5
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human LisH domain-containing protein ARMC9 protein (401-665AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:2000
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image



Description

Western blot All lanes: ARMC9 antibody at 6µg/ml Lane 1: Jurkat whole cell lysate Lane 2: CEM whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 92, 76 kDa Observed band size: 76 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO37430 at dilution of 1:100