

KATNAL2 Antibody

PACO37634

Description

This KATNAL2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO37634
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	DKFZP667C165 antibody, Katanin p60 ATPase-containing subunit A-like 2 antibody, Katanin p60 subunit A like 2 antibody, Katanin p60 subunit A-like 2 antibody, KATL2_HUMAN antibody, KATNAL 2 antibody, katnal2 antibody, MGC33211 antibody, p60 katanin like 2 antibody, p60 katanin-like 2 antibody
Clone:	Polyclonal
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q8IYT4
Host Species:	Rabbit
Purification:	Antigen Affinity Purified
Immunogen:	Recombinant Human Katanin p60 ATPase-containing subunit A-like 2 protein (18-307AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:2000
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image



Description

Western blot All lanes: KATNAL2 antibody at 1µg/ml + U251 whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 62, 53 kDa Observed band size: 62 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using PACO37634 at dilution of 1:100