

L3MBTL4 Antibody

PACO38562

Description

This L3MBTL4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO38562
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	L3MBTL4 antibody, Lethal(3)malignant brain tumor-like protein 4 antibody, H-l(3)mbt-like protein 4 antibody, L(3)mbt-like protein 4 antibody, L3mbt-like 4 antibody
Clone:	Polyclonal
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q8NA19
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Lethal(3)malignant brain tumor-like protein 4 protein (1-270AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data


Image



Description

Immunohistochemistry of paraffin-embedded human testis tissue using PACO38562 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human prostate cancer using PACO38562 at dilution of 1:100



Western Blot Positive WB detected in:
 Rat brain tissue All lanes: L3MBTL4
 antibody at 3µg/ml Secondary Goat
 polyclonal to rabbit IgG at 1/50000
 dilution Predicted band size: 72, 62 kDa
 Observed band size: 72 kDa