

ST6GALNAC3 Antibody

PACO38574

Description

This ST6GALNAC3 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO38574
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	ST6GALNAC3 antibody, SIAT7C antibody, UNQ2787/PRO7177Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 3 antibody, EC 2.4.99.7 antibody, GalNAc alpha-2,6-sialyltransferase III antibody, ST6GalNAc III antibody, ST6GalNAcIII antibody, STY antibody, Sialyltransferase 7C antibody, SIAT7-C antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q8NDV1
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 3 protein (40-305AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-