

UPRT Antibody

PACO38698

Description

This UPRT Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO38698
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	DKFZp781E1243 antibody, FUR1 antibody, Gm774 antibody, MGC23937 antibody, RGD1564806 antibody, RP11-311P8.3 antibody, RP23-92C4.1 antibody, UPP antibody, UPP_HUMAN antibody, uprt antibody, Uracil phosphoribosyltransferase homolog antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q96BW1
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Uracil phosphoribosyltransferase homolog protein (1-309AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human cervical cancer using PACO38698 at dilution of 1:100