

## CCDC127 Antibody

PACO38938

### Description

---

This CCDC127 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |   |
|----------------------|---|
| <b>SKU:</b>          | PACO38938   |
| <b>Contents:</b>     | 50µg<br>Bradford Reagent: 1 vial (2ml)                    |
| <b>Category:</b>     | -   |
| <b>Synonyms:</b>     | CCDC127Coiled-coil domain-containing protein 127 antibody |
| <b>Clone:</b>        | Polyclonal  |
| <b>Applications:</b> | <b>ELISA</b> <b>WB</b>                                    |
| <b>Conjugation:</b>  | Non-conjugated  |
| <b>Reactivity:</b>   | Human, Mouse  |

### Antibody Data

---

|                           |   |
|---------------------------|---|
| <b>Isotype:</b>           | IgG   |
| <b>Uniprot:</b>           | Q96BQ5  |
| <b>Host Species:</b>      | Rabbit  |
| <b>Purification:</b>      | >95%, Protein G purified  |
| <b>Immunogen:</b>         | Recombinant Human Coiled-coil domain-containing protein 127 protein (1-260AA) |
| <b>Immunogen Species:</b> | Homo sapiens (Human)  |
| <b>Buffer:</b>            | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| <b>Form:</b>              | Liquid  |

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.


| <b>Recommended Dilutions:</b> | <b>Application</b> | <b>Recommended Dilution</b> |
|-------------------------------|--------------------|-----------------------------|
|                               | WB                 | 1:500-1:2000                |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

### Image



### Description

Western blot All lanes: CCDC127 antibody at 6µg/ml + Mouse kidney tissue Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 31 kDa Observed band size: 31 kDa