

TSPAN14 Antibody

PACO39462

Description

This TSPAN14 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO39462
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	TSPAN14, TM4SF14, Tetraspanin-14, Tspan-14, DC-TM4F2, Transmembrane 4 superfamily member 14
Clone:	Polyclonal
Applications:	ELISA IHC IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q8NG11
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Tetraspanin-14 protein (114-215AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image




Description

Immunofluorescent analysis of U251 cells using PACO39462 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO39462 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39462 at dilution of 1:100