

sodB Antibody

PACO39770

Description

This sodB Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO39770
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	sodB antibody, AFUA_1G14550 antibody, Superoxide dismutase [Mn] antibody, mitochondrial antibody, EC 1.15.1.1 antibody, allergen Asp f 6 antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Neosartorya fumigata

Antibody Data

Isotype:	IgG
Uniprot:	Q92450
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Neosartorya fumigata Superoxide dismutase [Mn], mitochondrial protein (1-210AA)
Immunogen Species:	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-