

KRT35 Antibody

PACO40054

Description

This KRT35 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO40054
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Hair keratin antibody, Hair keratin, type I Ha5 antibody, HHA5 antibody, HKA5 antibody, K35 antibody, Keratin antibody, Keratin, type I cuticular Ha5 antibody, Keratin-35 antibody, KRT35 antibody, KRT35_HUMAN antibody, KRTHA5 antibody, type I cuticular Ha5 antibody, type I Ha5 antibody
Clone:	Polyclonal
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q92764
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Keratin, type I cuticular Ha5 protein (1-455AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image




Description

Western blot All lanes: KRT35 antibody at 4µg/ml + PC-3 whole cell lysate
Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 51 kDa Observed band size: 51 kDa



Immunohistochemistry of paraffin-embedded human skin tissue using PACO40054 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO40054 at dilution of 1:100