

CABP5 Antibody

PACO40626

Description

This CABP5 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO40626
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	CABP3 antibody, CaBP5 antibody, CABP5_HUMAN antibody, Calcium binding protein 3 antibody, Calcium-binding protein 5 antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q9NP86
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Calcium-binding protein 5 protein (1-173AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western Blot All lanes: CABP5 antibody at 12µg/ml Lane 1: Mouse kidney tissue Lane 2: Mouse thymus tissue Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 20 kDa Observed band size: 20 kDa