

PLEKHJ1 Antibody

PACO40886

Description

This PLEKHJ1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO40886
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	9530063M10Rik antibody, FLJ10297 antibody, Guanine nucleotide releasing protein x antibody, Guanine nucleotide-releasing protein x antibody, Likely ortholog of mouse guanine nucleotide releasing protein x antibody, PH domain containing family J member 1 antibody, PH domain-containing family J member 1 antibody, PKHJ1_HUMAN antibody, Pleckstrin homology domain containing family J member 1 antibody, Pleckstrin homology domain-containing family J member 1 antibody, PLEKHJ 1 antibody, PLEKHJ1 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9NW61
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Pleckstrin homology domain-containing family J member 1 protein (1-149AA)
Immunogen Species:	Homo sapiens (Human)

Buffer: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS,
PH 7.4

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human testis tissue using PACO40886 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human gastric cancer using PACO40886 at dilution of 1:100