

CDADC1 Antibody

PACO41266

Description

This CDADC1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO41266
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	bA103J18.1 antibody, CDAC1_HUMAN antibody, CDADC1 antibody, cytidine and dCMP deaminase domain containing 1 antibody, Cytidine and dCMP deaminase domain-containing protein 1 antibody, MGC150615 antibody, MGC41774 antibody, MGC57136 antibody, NYD-SP15 antibody, OTTHUMP00000018405 antibody, protein kinase NYD-SP15 antibody, RP11-103J18.2 antibody, Testis development protein NYD-SP15 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9BWV3
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Cytidine and dCMP deaminase domain-containing protein 1 protein (1-99AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO41266 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human testis tissue using PACO41266 at dilution of 1:100