

REEP4 Antibody

PACO41338

Description

This REEP4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO41338
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	REEP4, C8orf20, PP432, Receptor expression-enhancing protein 4
Clone:	Polyclonal
Applications:	ELISA WB IHC IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9H6H4
Host Species:	Rabbit
Purification:	Antigen Affinity Purified
Immunogen:	Recombinant Human Receptor expression-enhancing protein 4 protein (78-257AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000
	IHC	1:20-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image




Description


Western blot All lanes: REEP4 antibody at 0.7µg/ml + U251 whole cell lysate Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 30, 20 kDa Observed band size: 30 kDa




Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO41338 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using PACO41338 at dilution of 1:100



Immunofluorescent analysis of HeLa cells using PACO41338 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot Positive WB detected in: U251 whole cell lysate, A549 whole cell lysate, MCF-7 whole cell lysate, PC-3 whole cell lysate, 293 whole cell lysate
 All lanes: REEP4 antibody at 1:1500
 Secondary Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 30, 20 kDa
 Observed band size: 30 kDa