

VAT1L Antibody

PACO41362

Description

This VAT1L Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO41362
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Probable oxidoreductase KIAA1576 antibody, Synaptic vesicle membrane protein VAT 1 homolog like antibody, Synaptic vesicle membrane protein VAT-1 homolog-like antibody, VAT 1L antibody, VAT1L antibody, VAT1L_HUMAN antibody, Vesicle amine transport protein 1 homolog (T. californica) like antibody
Clone:	Polyclonal
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9HCJ6
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Synaptic vesicle membrane protein VAT-1 homolog-like protein (1-419AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:


Application	Recommended Dilution
WB	1:500-1:2000
IHC	1:20-1:200

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image



Description

Western blot All lanes: VAT1L antibody at 0.3µg/ml + U87 whole cell lysate
Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 46 kDa Observed band size: 46 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using PACO41362 at dilution of 1:100