

SEC14L4 Antibody

PACO41774

Description

This SEC14L4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO41774
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	AI256582 antibody, AW536553 antibody, dJ130H16 antibody, MGC30313 antibody, RGD1565810 antibody, S14L4_HUMAN antibody, SEC14, S. cerevisiae, homolog of, 4 antibody, SEC14-like 4 (S. cerevisiae) antibody, SEC14-like 4 antibody, SEC14-like protein 4 antibody, Sec14l4 antibody, SEC14p-like protein TAP3 antibody, SPF antibody, TAP antibody, TAP3 antibody, Tocopherol-associated protein 3 antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9UDX3
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human SEC14-like protein 4 protein (201-360AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human kidney tissue using PACO41774 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human liver tissue using PACO41774 at dilution of 1:100