

ATP6V0A2 Antibody

PACO42834

Description

This ATP6V0A2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO42834
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	a2 antibody, A2V ATPase antibody, ARCL antibody, ATP6a2 antibody, ATP6N1D antibody, ATP6V0A2 antibody, ATPase, H ⁺ transporting, lysosomal V0 subunit a isoform 2 antibody, ATPase, H ⁺ transporting, lysosomal V0 subunit a2 antibody, Infantile malignant osteopetrosis antibody, J6B7 antibody, Lysosomal H(+) transporting ATPase V0 subunit a2 antibody, Lysosomal H(+)-transporting ATPase V0 subunit a2 antibody, regeneration and tolerance factor antibody, Stv1 antibody, TJ6 antibody, TJ6M antibody, TJ6s antibody, V ATPase 116 kDa isoform a2 antibody, V type proton ATPase 116 kDa subunit a antibody, V t...Read more
Clone:	Polyclonal
Applications:	ELISA IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q9Y487
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human V-type proton ATPase 116 kDa subunit a isoform 2 protein (476-549AA)
Immunogen Species:	Homo sapiens (Human)

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Buffer: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunofluorescent analysis of A431 cells using PACO42834 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)