

ANKS6 Antibody

PACO46594

Description

This ANKS6 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO46594
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	ANKRD14 antibody, Anks6 antibody, ANKS6_HUMAN antibody, ankyrin repeat and SAM domain containing protein 6 antibody, Ankyrin repeat and SAM domain-containing protein 6 antibody, Ankyrin repeat and sterile alpha motif domain containing 6 antibody, ankyrin repeat domain 14 antibody, ankyrin repeat domain containing protein 14 antibody, Ankyrin repeat domain-containing protein 14 antibody, SAM domain containing protein 6 antibody, SAM domain-containing protein 6 antibody, SamCystin antibody, SAMD6 antibody, sterile alpha motif domain containing 6 antibody, sterile alpha motif domain containing pr... Read more
Clone:	Polyclonal
Applications:	ELISA IHC IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q68DC2
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Ankyrin repeat and SAM domain-containing protein 6 protein (1-132AA)
Immunogen Species:	Homo sapiens (Human)

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Buffer: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS,
PH 7.4

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO46594 at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using PACO46594 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)