

DCAF8 Antibody

PACO46602

Description

This DCAF8 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO46602
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	DCAF8 antibody, H326 antibody, WDR42ADDB1- and CUL4-associated factor 8 antibody, WD repeat-containing protein 42A antibody
Clone:	Polyclonal
Applications:	ELISA WB IHC
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q5TAQ9
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human DDB1- and CUL4-associated factor 8 protein (1-135AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:2000-1:5000
	IHC	1:20-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image




Description

Western Blot Positive WB detected in: Mouse heart tissue, Mouse lung tissue
All lanes: DCAF8 antibody at 3.5µg/ml
Secondary Goat polyclonal to rabbit IgG at 1/50000 dilution Predicted band size: 67, 31 kDa Observed band size: 67 kDa



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO46602 at dilution of 1:100



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO46602 at dilution of 1:100