

GALNT11 Antibody

PACO46946

Description

This GALNT11 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO46946
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	A430075I06Rik antibody, A1648252 antibody, E430002F06Rik antibody, FLJ21634 antibody, GalNAc Transferase 11 antibody, GalNAc-T11 antibody, GalNAcT11 antibody, GALNT11 antibody, GLT11_HUMAN antibody, MGC71630 antibody, Polypeptide GalNAc transferase 11 antibody, Polypeptide N-acetylgalactosaminyltransferase 11 antibody, pp-GaNTase 11 antibody, Protein-UDP acetylgalactosaminyltransferase 11 antibody, tcag7.1057 antibody, UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 11 antibody, UDP-N-acetyl-alpha-D-galactosamine:polypeptide N-acetylgalactosaminyltransferase 11 (GalNAc-T11) antibody
Clone:	Polyclonal
Applications:	ELISA IHC IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q8NCW6
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 11 protein (420-606AA)

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Immunogen Species: Homo sapiens (Human)

Buffer: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data


Image



Description

Immunohistochemistry of paraffin-embedded human gastric cancer using PACO46946 at dilution of 1:100

Immunohistochemistry of paraffin-embedded human cervical cancer using PACO46946 at dilution of 1:100



Immunofluorescent analysis of HeLa cells using PACO46946 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)