

FBXL21 Antibody

PACO47230

Description

This FBXL21 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO47230
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	F box and leucine rich repeat protein 21 antibody, F box and leucine rich repeat protein 3 pseudogene antibody, F box and leucine rich repeat protein 3B antibody, F box protein Fbl3b antibody, F box/LRR repeat protein 21 antibody, F box/LRR repeat protein 3B antibody, F-box and leucine-rich repeat protein 21 antibody, F-box and leucine-rich repeat protein 3B antibody, F-box/LRR-repeat protein 21 antibody, F-box/LRR-repeat protein 3B antibody, Fbl 21 antibody, FBL 3 antibody, FBL 3B antibody, Fbl21 antibody, FBL3 antibody, FBL3B antibody, FBXL 21 antibody, FBXL 3B antibody, FBXL 3P antibody, Fb... Read more
Clone:	Polyclonal
Applications:	ELISA WB IHC IF
Conjugation:	Non-conjugated
Reactivity:	Human, Mouse

Antibody Data

Isotype:	IgG
Uniprot:	Q9UKT6
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human F-box/LRR-repeat protein 21 protein (78-137AA)
Immunogen Species:	Homo sapiens (Human)

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Buffer: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:1000-1:5000
	IHC	1:20-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol


Validation Data

Image




Description

Western blot All lanes: FBL21 antibody at 10 µg/ml
Lane 1: Mouse heart tissue Lane 2: Mouse skeletal muscle tissue Secondary Goat polyclonal to rabbit IgG at 1/10000 dilution Predicted band size: 49 kDa Observed band size: 49 kDa



Immunohistochemistry of paraffin-embedded human lung cancer using PACO47230 at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using PACO47230 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)