

GRHL3 Antibody

PACO49258

Description

This GRHL3 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO49258
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Drosophila antibody, Grainyhead like 3 antibody, Grainyhead like protein 3 homolog antibody, Grainyhead-like protein 3 homolog antibody, grhl3 antibody, GRHL3_HUMAN antibody, Sister of mammalian grainyhead antibody, Sister of mammalian grainyhead protein antibody, SOM antibody, TFCP2L4 antibody, Transcription factor CP2 like 4 antibody, Transcription factor CP2-like 4 antibody
Clone:	Polyclonal
Applications:	ELISA IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q8TE85
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Grainyhead-like protein 3 homolog protein (31-182AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunofluorescent analysis of PC-3 cells using PACO49258 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)