

## CDH-1 Antibody

PACO50254

### Description

---

This CDH-1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	PACO50254
<b>Contents:</b>	50µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	-
<b>Synonyms:</b>	CDH-1, antibody, CDH-2Cellobiose dehydrogenase antibody, CDH antibody, EC 1.1.99.18 antibody, Cellobiose-quinone oxidoreductase antibody
<b>Clone:</b>	Polyclonal
<b>Applications:</b>	<a href="#">ELISA</a>
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	Phanerochaete chrysosporium

### Antibody Data

---

<b>Isotype:</b>	IgG
<b>Uniprot:</b>	Q01738
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	>95%, Protein G purified
<b>Immunogen:</b>	Recombinant Phanerochaete chrysosporium Cellobiose dehydrogenase protein (19-208AA)
<b>Immunogen Species:</b>	Phanerochaete chrysosporium (White-rot fungus) (Sporotrichum pruinosum)
<b>Buffer:</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Form:</b>	Liquid

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

Image	Description
-	-