

ROS1 Antibody

PACO51066

Description

This ROS1 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO51066
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	ROS1 antibody, DML1 antibody, At2g36490 antibody, F1O11.12 antibody, DNA glycosylase/AP lyase ROS1 antibody, EC 4.2.99.18 antibody, DEMETER-like protein 1 antibody, Protein REPRESSOR OF SILENCING 1 antibody, Protein ROS1 antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Arabidopsis thaliana

Antibody Data

Isotype:	IgG
Uniprot:	Q9SJQ6
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Arabidopsis thaliana Protein ROS1 protein (1-320AA)
Immunogen Species:	Arabidopsis thaliana (Mouse-ear cress)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-