

Tropomyosin Antibody

PACO51154

Description

This Tropomyosin Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO51154
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Tropomyosin antibody, Allergen Cha f I antibody, allergen Cha f 1 antibody, Fragment antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Charybdis feriata

Antibody Data

Isotype:	IgG
Uniprot:	Q9N2R3
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Charybdis feriata Tropomyosin protein (1-264AA)
Immunogen Species:	Charybdis feriata (Crucifix crab) (Cancer feriatus)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-