

mrcA Antibody

PACO52254

Description

This mrcA Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO52254
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	mrcA antibody, ponA antibody, NMB1807 antibody, Penicillin-binding protein 1A antibody, PBP-1a antibody, PBP1a) [Includes: Penicillin-insensitive transglycosylase antibody, EC 2.4.1.129 antibody, Peptidoglycan TGase), Penicillin-sensitive transpeptidase antibody, EC 3.4.16.4 antibody, DD-transpeptidase]] antibody
Clone:	Polyclonal
Applications:	ELISA
Conjugation:	Non-conjugated
Reactivity:	Neisseria meningitidis serogroup B

Antibody Data

Isotype:	IgG
Uniprot:	POA0Z6
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Neisseria meningitidis serogroup B Penicillin-binding protein 1A protein (206-413AA)
Immunogen Species:	Neisseria meningitidis serogroup B (strain MC58)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage

Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	-	-

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image	Description
-	-