

Loxl4 Antibody

PACO52658

Description

This Loxl4 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO52658
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Loxl4 antibody, Loxc antibody, Lysyl oxidase homolog 4 antibody, EC 1.4.3.-antibody, Lysyl oxidase-like protein 4 antibody, Lysyl oxidase-related protein C antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Mouse, Rat

Antibody Data

Isotype:	IgG
Uniprot:	Q924C6
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Mouse Lysyl oxidase homolog 4 protein (26-757AA)
Immunogen Species:	Mus musculus (Mouse)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western Blot Positive WB detected in: Rat brain tissue All lanes: Loxl4 antibody at 4µg/ml Secondary Goat polyclonal to rabbit IgG at 1/50000 dilution Predicted band size: 85 kDa Observed band size: 85 kDa