

TFDP3 Antibody

PACO53474

Description

This TFDP3 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO53474
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Cancer/testis antigen 30 antibody, CT30 antibody, DP4 antibody, E2F like protein antibody, E2F-like antibody, HCA661 antibody, Hepatocellular carcinoma-associated antigen 661 antibody, TFDP3 antibody, TFDP3_HUMAN antibody, Transcription factor Dp family member 3 antibody
Clone:	Polyclonal
Applications:	ELISA IHC IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	Q5H9I0
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human Transcription factor Dp family member 3 protein (61-196AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:20-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunohistochemistry of paraffin-embedded human testis tissue using PACO53474 at dilution of 1:100

Immunofluorescent analysis of HepG2 cells using PACO53474 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)