

Pcsk2 Antibody

PACO53710

Description

This Pcsk2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO53710
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Pcsk2 antibody, Nec-2 antibody, Nec2 antibody, Rpc2 antibody, Neuroendocrine convertase 2 antibody, NEC 2 antibody, EC 3.4.21.94 antibody, KEX2-like endoprotease 2 antibody, Prohormone convertase 2 antibody, Proprotein convertase 2 antibody, PC2 antibody
Clone:	Polyclonal
Applications:	ELISA WB
Conjugation:	Non-conjugated
Reactivity:	Rat, Human

Antibody Data

Isotype:	IgG
Uniprot:	P28841
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Rat Neuroendocrine convertase 2 protein (109-637AA)
Immunogen Species:	Rattus norvegicus (Rat)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	WB	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Western Blot Positive WB detected in: Jurkat whole cell lysate All lanes: Pcsk1 antibody at 3µg/ml Secondary Goat polyclonal to rabbit IgG at 1/50000 dilution Predicted band size: 71 kDa Observed band size: 71 kDa