

## notch1a Antibody

PACO54242

### Description

---

This notch1a Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	PACO54242
<b>Contents:</b>	50µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	-
<b>Synonyms:</b>	notch1a antibody, notch antibody, Neurogenic locus notch homolog protein 1 antibody, Notch 1) [Cleaved into: Notch 1 extracellular truncation antibody, NEXT), Notch 1 intracellular domain antibody, NICD)] antibody
<b>Clone:</b>	Polyclonal
<b>Applications:</b>	<b>ELISA</b>
<b>Conjugation:</b>	Non-conjugated
<b>Reactivity:</b>	Zebrafish

### Antibody Data

---

<b>Isotype:</b>	IgG
<b>Uniprot:</b>	P46530
<b>Host Species:</b>	Rabbit
<b>Purification:</b>	>95%, Protein G purified
<b>Immunogen:</b>	Recombinant Zebrafish Neurogenic locus notch homolog protein 1 protein (21-338AA)
<b>Immunogen Species:</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Buffer:</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Form:</b>	Liquid

## Preparation & Storage

---

**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

<b>Recommended Dilutions:</b>	<b>Application</b>	<b>Recommended Dilution</b>
	-	-

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

---

<b>Image</b>	<b>Description</b>
-	-