

## SAA2 Antibody

PACO58885

### Description

---

This SAA2 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |  |
|----------------------|--|
| <b>SKU:</b>          | PACO58885  |
| <b>Contents:</b>     | 50µg<br>Bradford Reagent: 1 vial (2ml)                           |
| <b>Category:</b>     | -  |
| <b>Synonyms:</b>     | SAA2 antibody, Serum amyloid A-2 protein antibody, SAA2 antibody |
| <b>Clone:</b>        | Polyclonal   |
| <b>Applications:</b> | <b>ELISA</b> <b>IHC</b>  |
| <b>Conjugation:</b>  | Non-conjugated   |
| <b>Reactivity:</b>   | Human  |

### Antibody Data

---

|                           |   |
|---------------------------|---|
| <b>Isotype:</b>           | IgG   |
| <b>Uniprot:</b>           | P0DJ19  |
| <b>Host Species:</b>      | Rabbit  |
| <b>Purification:</b>      | >95%, Protein G purified  |
| <b>Immunogen:</b>         | Recombinant Human Serum amyloid A-2 protein (19-94AA)                         |
| <b>Immunogen Species:</b> | Homo sapiens (Human)  |
| <b>Buffer:</b>            | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| <b>Form:</b>              | Liquid  |

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

## Preparation & Storage



**Storage:** Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  
Store Bradford Reagent at Room Temperature for 1 Year.

| Recommended Dilutions: | Application | Recommended Dilution |
|------------------------|-------------|----------------------|
|                        | IHC         | 1:500-1:1000         |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

## Validation Data

### Image



### Description

IHC image of PACO58885 diluted at 1:600 and staining in paraffin-embedded human liver cancer performed on a Leica Bond<sup>TM</sup> system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.

IHC image of PACO58885 diluted at 1:600 and staining in paraffin-embedded human liver tissue performed on a Leica Bond<sup>TM</sup> system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.