

TYW5 Antibody

PACO59684

Description

This TYW5 Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO59684
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	TYW5 antibody, C2orf60 antibody, tRNA wybutosine-synthesizing protein 5 antibody, hTYW5 antibody, EC 1.14.11.42 antibody, tRNA(Phe antibody, 7-(3-amino-3-carboxypropyl)wyosine(37)-C(2))-hydroxylase antibody
Clone:	Polyclonal
Applications:	ELISA IF
Conjugation:	Non-conjugated
Reactivity:	Human

Antibody Data

Isotype:	IgG
Uniprot:	A2RUC4
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Human tRNA wybutosine-synthesizing protein 5 protein (1-152AA)
Immunogen Species:	Homo sapiens (Human)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

Immunofluorescence staining of HeLa cells with PACO59684 at 1:166, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).