

Ace Antibody

PACO62663

Description

This Ace Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	PACO62663
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Category:	-
Synonyms:	Ace antibody, Dcp1 antibody, Angiotensin-converting enzyme antibody, ACE antibody, EC 3.2.1.- antibody, EC 3.4.15.1 antibody, Dipeptidyl carboxypeptidase I antibody, Kininase II antibody, CD antigen CD143) [Cleaved into: Angiotensin-converting enzyme antibody, soluble form] antibody
Clone:	Polyclonal
Applications:	ELISA IHC
Conjugation:	Non-conjugated
Reactivity:	Rat

Antibody Data

Isotype:	IgG
Uniprot:	P47820
Host Species:	Rabbit
Purification:	>95%, Protein G purified
Immunogen:	Recombinant Rat Angiotensin-converting enzyme protein (1046-1263AA)
Immunogen Species:	Rattus norvegicus (Rat)
Buffer:	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form:	Liquid

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Preparation & Storage


Storage: Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions:	Application	Recommended Dilution
	IHC	1:200-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Validation Data

Image



Description

IHC image of PACO62663 diluted at 1:400 and staining in paraffin-embedded rat kidney tissue performed on a Leica BondTM system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a biotinylated secondary antibody and visualized using an HRP conjugated SP system.