



EGTT-Scavenger for GSH/GSSG Assays (BA0246)

Product Datasheet | SKU BA0246

Description

A ready-to-use chemical reagent (500 µL, sufficient for 100 tests) used for sample treatment to scavenge reduced glutathione (GSH) from samples, enabling selective measurement of oxidised glutathione (GSSG) with the companion GSH/GSSG Assay Kit. Necessary for cell lysate samples.

Specifications

Species Reactivity	All
Kit Size	100 Assays
Equipment Required	Microplate reader
Storage	Room Temperature
Shelf Life	12 months
Shipping	Room Temperature

Kit Components

Component	Quantity	Storage
EGTT-Scavenger solution	500 µL	Room Temperature

Ordering Information

Assay Genie Cat. No.	Pack Size
BA0246	100 Assays

Protocol Summary

Prepare and use the reagent as described in the kit insert. To be used in conjunction with the companion GSH/GSSG (glutathione) Assay Kit and its associated materials.

Notes

- This reagent is an accessory intended for use with the companion GSH/GSSG (glutathione) Assay Kit.
- Reagents are for research use only. Refer to the MSDS for detailed safety information.

Assay Genie | sales@assaygenie.com | www.assaygenie.com

Technical specifications are based on manufacturer documentation and may be updated. Protocols shown are summaries - always follow the instructions supplied with your kit/lot.