

## Anti-Denosumab Neutralizing Antibody ADA ELISA Kit

SKU: GXEL00013 | Assay: Quantitative ELISA (Colorimetric)

The Anti-Denosumab Neutralizing Antibody ADA ELISA Kit (SKU: GXEL00013) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-Denosumab Neutralizing Antibody in Plasma, Serum. The assay covers a detection range of 15.625-1,000 ng/ml and is read colorimetrically at 450 nm.

### Specifications

<b>Product Name</b>	Anti-Denosumab Neutralizing Antibody ADA ELISA Kit
<b>SKU</b>	GXEL00013
<b>Accession</b>	O14788
<b>Alternative Names</b>	AMG162, HL2001, 615258-40-7
<b>Assay Type</b>	Quantitative
<b>Detection Method</b>	Colorimetric
<b>Detection Range</b>	15.625-1,000 ng/ml
<b>Recovery</b>	80-120%
<b>Sample Types</b>	Plasma, Serum
<b>Storage</b>	2-8°C. See Storage & Stability.
<b>Note</b>	For Research Use Only. Not for diagnostic procedures.

### Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00013.pdf> | Version 1.0, 09 July 2026