

Serclutamab ELISA Kit

SKU: GXEL00025 | Assay: Quantitative ELISA (Colorimetric)

The Serclutamab ELISA Kit (SKU: GXEL00025) is a high-sensitivity quantitative ELISA for the quantitative measurement of Serclutamab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Serclutamab talirine (ABBV-321) is an epidermal growth factor receptor (EGFR)-targeted antibody-drug conjugate that consists of a humanized immunoglobulin G1 anti-EGFR monoclonal antibody conjugated to a pyrrolobenzodiazepine (PBD) dimer via a maleimidocaproyl-valine-alanine linker.

Specifications

Product Name	Serclutamab ELISA Kit
SKU	GXEL00025
Accession	P00533
Alternative Names	ABBV-321, DC-1630423, PR-1594407, Serclutamab talirine, 2140172-41-2
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00025.pdf> | Version 1.0, 09 July 2026