

Zanolimumab ELISA Kit

SKU: GXEL00056 | Assay: Quantitative ELISA (Colorimetric)

The Zanolimumab ELISA Kit (SKU: GXEL00056) is a high-sensitivity quantitative ELISA for the quantitative measurement of Zanolimumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Zanolimumab is a biosimilar that targets CD4. The CD4 antigen is involved in the recognition of MHC class II molecules and is a co-receptor for HIV. CD4 is primarily expressed in a subset of T-lymphocytes, also referred to as T helper cells, but may also be expressed by other cells in the immune system, such as monocytes, macrophages, and dendritic cells. At the tissue level, CD4 expression may be detected in thymus, lymph nodes, tonsils, and spleen, and also in specific regions of the brain, gut, and other non-lymphoid tissues. CD4 functions to initiate or augment the early phase of T-cell activation through its association with the T-cell receptor complex and protein tyrosine kinase, Lck. It may also function as an important mediator of direct neuronal damage in infectious and immune-mediated diseases of the central nervous system. Multiple alternatively spliced transcripts have been identified in this gene [RefSeq, July 2017].

Specifications

Product Name	Zanolimumab ELISA Kit
SKU	GXEL00056
Accession	P01730
Alternative Names	HuMax-CD4, 652153-01-0
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00056.pdf> | Version 1.0, 09 July 2026