

Palivizumab ELISA Kit

SKU: GXEL00062 | Assay: Quantitative ELISA (Colorimetric)

The Palivizumab ELISA Kit (SKU: GXEL00062) is a high-sensitivity quantitative ELISA for the quantitative measurement of Palivizumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Palivizumab (brand name Synagis) is a humanized IgG monoclonal antibody (mAb) manufactured by MedImmune with recombinant DNA technology. It is used in the prevention of respiratory syncytial virus (RSV) infections. It is recommended for infants that are high-risk because of prematurity or other medical problems such as congenital heart disease. Palivizumab directed against an epitope in the A antigenic site of the F protein of RSV. Palivizumab is a composite of human (95%) and murine (5%) antibody sequences. The human heavy chain sequence was derived from the constant domains of human IgG1 and the variable framework regions of the VH genes Cor and Cess. The human light chain sequence was derived from the constant domain of Cκ and the variable framework regions of the VL gene K104 with Jκ -4. The murine sequences were derived from a murine monoclonal antibody, Mab 1129, in a process that involved the grafting of the murine complementarity determining regions into the human antibody frameworks. Palivizumab is composed of two heavy chains and two light chains and has a molecular weight of approximately 148, 000 Daltons. Palivizumab was approved for medical use in 1998 but the indications for prophylaxis vary worldwide. Palivizumab is primarily indicated for the following conditions: prematurity (gestational age ≤ 35 weeks), bronchopulmonary dysplasia/chronic lung disease, hemodynamically significant congenital heart disease, and other serious medical disorders on an ad hoc basis. Infants with these medical conditions are at increased risk of more severe disease and more serious subsequent sequelae.

Specifications

Product Name	Palivizumab ELISA Kit
SKU	GXEL00062
Accession	P03420
Alternative Names	MEDI-493, 188039-54-5, MEDI493, SYNAGIS, Humanized MAb 1129, 1129
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00062.pdf> | Version 1.0, 09 July 2026