

Namulumab ELISA Kit

SKU: GXEL00067 | Assay: Quantitative ELISA (Colorimetric)

The Namulumab ELISA Kit (SKU: GXEL00067) is a high-sensitivity quantitative ELISA for the quantitative measurement of Namulumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Namulumab (alternative identifier MT203) is a human monoclonal antibody (class IgG1 kappa) that targets granulocyte macrophage-colony stimulating factor (GM-CSF), which is also known as colony stimulating factor 2 (CSF2). Namulumab was produced by Micromet Inc and is under development by Takeda Pharmaceuticals International. It is currently being researched for application in rheumatoid arthritis (RA) and psoriatic arthritis. Clinical trials investigating the therapeutic utility of Namulumab have include phase I and phase II clinical trials to establish the safety, tolerability and preliminary therapeutic utility of the antibody in plaque psoriasis and rheumatoid arthritis.

Specifications

Product Name	Namulumab ELISA Kit
SKU	GXEL00067
Accession	P04141
Alternative Names	MT203, 1206681-39-1
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00067.pdf> | Version 1.0, 09 July 2026