

Zanidatamab ELISA Kit

SKU: GXEL00073 | Assay: Quantitative ELISA (Colorimetric)

The Zanidatamab ELISA Kit (SKU: GXEL00073) is a high-sensitivity quantitative ELISA for the quantitative measurement of Zanidatamab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: zanidatamab is a novel her2-targeted antibody that binds two distinct extracellular domains of her2, allowing for multiple mechanisms of action including enhanced binding, clustering, receptor internalization and downregulation

Specifications

Product Name	Zanidatamab ELISA Kit
SKU	GXEL00073
Accession	P04626
Alternative Names	ZW-25, ZW25, zanidatamab zovodotin, 2169946-15-8
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00073.pdf> | Version 1.0, 09 July 2026