

Sirukumab ELISA Kit

SKU: GXEL00080 | Assay: Quantitative ELISA (Colorimetric)

The Sirukumab ELISA Kit (SKU: GXEL00080) is a high-sensitivity quantitative ELISA for the quantitative measurement of Sirukumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Sirukumab (developmental code name CNTO-136) is a human monoclonal antibody (mAb) designed for the treatment of rheumatoid arthritis (RA). It acts against the proinflammatory cytokine interleukin (IL) -6 (IL-6). Sirukumab was evaluated in five Phase 3 studies of patients with RA. In September 2016, Janssen submitted a biologics license application (BLA) seeking approval of sirukumab for the treatment of moderately to severely active RA, but they announced in September 2017 that they received a complete response letter from the Food and Drug Administration (FDA) indicating additional clinical data were needed to further evaluate the safety of sirukumab in that patient population. Janssen subsequently decided not to pursue global approvals of sirukumab for the treatment of moderately to severely active RA. Sirukumab was being evaluated in a Phase 3 study (NCT02531633) of another inflammatory disease, giant cell arteritis. This study, sponsored by GlaxoSmithKline (GSK), was terminated in November 2017 due to GSK's decision to return rights to sirukumab to Janssen and discontinue sirukumab development in giant cell arteritis. Also, in November 2017, a Phase 3 study (NCT02899026) sponsored by GSK of sirukumab in patients with polymyalgia rheumatica, an inflammatory disorder that causes muscle pain and stiffness, was withdrawn prior to enrollment.

Specifications

Product Name	Sirukumab ELISA Kit
SKU	GXEL00080
Accession	P05231
Alternative Names	CNTO 136, 1194585-53-9
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00080.pdf> | Version 1.0, 09 July 2026