

Anti-Bevacizumab Neutralizing Antibody ADA ELISA Kit

SKU: GXEL00135 | Assay: Quantitative ELISA (Colorimetric)

The Anti-Bevacizumab Neutralizing Antibody ADA ELISA Kit (SKU: GXEL00135) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-Bevacizumab Neutralizing Antibody in Plasma, Serum. The assay covers a detection range of 15.625-1, 000 ng/ml and is read colorimetrically at 450 nm.

Specifications

Product Name	Anti-Bevacizumab Neutralizing Antibody ADA ELISA Kit
SKU	GXEL00135
Accession	P15692
Alternative Names	12-IgG1, F(ab)-12 IgG1, Fab-12 IgG1, rhuMAb-VEGF, ABP 215, 216974-75-3
Assay Type	Quantitative
Detection Method	Colorimetric
Detection Range	15.625-1, 000 ng/ml
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00135.pdf> | Version 1.0, 09 July 2026