

## Anti-Daratumumab ADA ELISA Kit

SKU: GXEL00171 | Assay: Quantitative ELISA (Colorimetric)

The Anti-Daratumumab ADA ELISA Kit (SKU: GXEL00171) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-Daratumumab in Plasma, Serum. The assay covers a detection range of 156.25 - 10, 000 ng/mL and is read colorimetrically at 450 nm.

### Specifications

|                          |                                                       |
|--------------------------|-------------------------------------------------------|
| <b>Product Name</b>      | Anti-Daratumumab ADA ELISA Kit                        |
| <b>SKU</b>               | GXEL00171                                             |
| <b>Accession</b>         | P28907                                                |
| <b>Alternative Names</b> | 945721-28-8                                           |
| <b>Assay Type</b>        | Quantitative                                          |
| <b>Detection Method</b>  | Colorimetric                                          |
| <b>Detection Range</b>   | 156.25 - 10, 000 ng/mL                                |
| <b>Recovery</b>          | 80-120%                                               |
| <b>Sample Types</b>      | Plasma, Serum                                         |
| <b>Storage</b>           | 2-8°C. See Storage & Stability.                       |
| <b>Note</b>              | For Research Use Only. Not for diagnostic procedures. |

### Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00171.pdf> | Version 1.0, 09 July 2026