

Ustekinumab ELISA Kit

SKU: GXEL00175 | Assay: Quantitative ELISA (Colorimetric)

The Ustekinumab ELISA Kit (SKU: GXEL00175) is a high-sensitivity quantitative ELISA for the quantitative measurement of Ustekinumab in Plasma, Serum. The assay covers a detection range of 31.25 - 2, 000 ng/mL and is read colorimetrically at 450 nm.

Specifications

Product Name	Ustekinumab ELISA Kit
SKU	GXEL00175
Accession	P29460
Alternative Names	CNTO 1275, TT-20, 815610-63-0
Assay Type	Quantitative
Detection Method	Colorimetric
Detection Range	31.25 - 2, 000 ng/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00175.pdf> | Version 1.0, 09 July 2026