

Cusatuzumab ELISA Kit

SKU: GXEL00178 | Assay: Quantitative ELISA (Colorimetric)

The Cusatuzumab ELISA Kit (SKU: GXEL00178) is a high-sensitivity quantitative ELISA for the quantitative measurement of Cusatuzumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Cusatuzumab is designed as a first-in-class, high-affinity anti-CD70 monoclonal antibody with multiple intended mechanisms of action, including Fc-mediated cytotoxicity with enhanced ADCC and inhibition of CD70/CD27 signaling, resulting in killing of malignant cells such as leukemia blasts and stem cells.

Specifications

Product Name	Cusatuzumab ELISA Kit
SKU	GXEL00178
Accession	P32970
Alternative Names	ARGX-110, 1399546-01-0
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00178.pdf> | Version 1.0, 09 July 2026