

Anetumab ELISA Kit

SKU: GXEL00207 | Assay: Quantitative ELISA (Colorimetric)

The Anetumab ELISA Kit (SKU: GXEL00207) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anetumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Anetumab ravtansine, also known as BAY 94-9343, is a novel, selective and highly potent antibody-drug conjugate consisting of a human anti-mesothelin antibody conjugated to the maytansinoid tubulin inhibitor DM4 via a disulfide-containing linker (a reducible SPDB linker [N-succinimidyl 4- (2-pyridyldithio)butanoate]) on an average of 3 lysyl.

Specifications

Product Name	Anetumab ELISA Kit
SKU	GXEL00207
Accession	Q13421
Alternative Names	BAY 94-9343, BAY 86-1903, BAY 2287409, Anetumab Ravtansine, 1954758-84-9
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00207.pdf> | Version 1.0, 09 July 2026