

Anti-Pembrolizumab Neutralizing Antibody ADA ELISA Kit

SKU: GXEL00217 | Assay: Quantitative ELISA (Colorimetric)

The Anti-Pembrolizumab Neutralizing Antibody ADA ELISA Kit (SKU: GXEL00217) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-Pembrolizumab Neutralizing Antibody in Plasma, Serum. The assay covers a detection range of 156.25 - 10, 000 ng/mL with a sensitivity of 83.4 ng/mL. and is read colorimetrically at 450 nm.

Specifications

Product Name	Anti-Pembrolizumab Neutralizing Antibody ADA ELISA Kit
SKU	GXEL00217
Accession	Q15116
Alternative Names	MK-3475
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	83.4 ng/mL.
Detection Range	156.25 - 10, 000 ng/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at -20°C, the rest reagents should be store at 4°C. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://store-h68l9z2lhx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00217.pdf> | Version 1.0, 09 July 2026