

Ixekizumab ELISA Kit

SKU: GXEL00221 | Assay: Quantitative ELISA (Colorimetric)

The Ixekizumab ELISA Kit (SKU: GXEL00221) is a high-sensitivity quantitative ELISA for the quantitative measurement of Ixekizumab in Plasma, Serum. The assay covers a detection range of 78.13 - 5, 000 ng/mL with a sensitivity of 39.3 ng/mL and is read colorimetrically at 450 nm.

Background: Ixekizumab is a humanized monoclonal antibody to interleukin-17A which acts as an antiinflammatory agent and is used to treat moderate-to-severe plaque psoriasis. Ixekizumab has not been linked to serum enzyme elevations during therapy or to instances of idiosyncratic acute liver injury.

Specifications

Product Name	Ixekizumab ELISA Kit
SKU	GXEL00221
Accession	Q16552
Alternative Names	LY2439821, 1143503-69-8
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	39.3 ng/mL
Detection Range	78.13 - 5, 000 ng/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at -20°C, the rest reagents should be store at 4°C. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00221.pdf> | Version 1.0, 09 July 2026